



# POOLED COMPOSITE & SINGLE-SITE EXAMPLE SHEET

## SURFACE

**1 Single-Site  
Surface Sample**

=



5 Surfaces in one location with 1 swab.

**Example:** Office space 1: Desk, phone, keyboard, Computer mouse, and light switch.

**1 Pooled Surface  
Sample**

=



Combine up to 5 Single-Site Samples  
(Total of 25 surfaces).

## AIR

**1 Single-Site  
Air Sample**

=



Up to 5 air sources in one location with 1 cassette.

**Example:** Front lobby: near entrance, welcome desk, bathroom door, vending machine, and water fountain.

**1 Pooled  
Air Sample**

=



Combine up to 5 Single-Site Samples  
(Total of 25 air sources).

**NOTE:** Both **air** & **surface** samples  
can be pooled together.